




neodisher® SBK



Special rinse aid for use in washer-disinfectors for bedpans



Liquid concentrate

Fields of application:

- Water softening and rinsing when reprocessing bedpans and urine bottles in bedpan washer-disinfectors with steam disinfection

Performance spectrum:

- Effectively prevents lime residues in steam generators, pipes and jets also at high degrees of water hardness
- With special softening components for effectively complexing water hardness
- For rapid and stain-free drying

Special properties:

- Foam-reducing; suitable for all acid-compatible washer-disinfectors for bedpans and urine bottles with steam disinfection
- Suitable for all degrees of water hardness
- Acidic, based on organic acids

Application and dosage:

neodisher SBK is used in bedpan washer-disinfectors with steam disinfection. neodisher SBK is dosed using dosing equipment installed in the washer-disinfector. The dosing amount is 0.5 - 1.5 ml/l depending on water hardness. Dosing systems and steam generators must be acid-compatible.

For removing stubborn soiling we recommend the additional use of neodisher SBR extra, a special detergent for bedpans and urine bottles.

Notes on application:

- For professional use only.
- Rinse out dosing system including suction hose with water before changing product.

- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.
- Please observe the recommendations given by the manufacturer of the items to be reprocessed.
- The instructions given by the manufacturer of the washer disinfector are to be observed.
- Do not mix with other products.

Technical data:

pH-range	3.7 - 3.0 (0.5 - 1.5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	1.1 g/cm ³ (20 °C)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

< 5 % polycarboxylates
also preservatives (potassium sorbate, methylchloroisothiazolinone/methylisothiazolinone)

CE-marks:

neodisher SBK conforms to the European directive 93/42/EEC, Annex I, concerning medical devices.

Storage information:

Always store at a temperature between 0 and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol





neodisher[®] SBK

Hazard and precautionary statements:

For safety information see EC safety data sheets.
These are available e.g. at www.drweigert.com under
the headline "Service".

Dispose only when container is empty and closed.
For disposal of product residues, refer to the
Material Safety Data Sheet.

MB 4006/3-2 Date of issue: 01/2016



With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.
