

Ordering Guide








A large, stylized graphic of a blue wave or series of concentric arcs, rendered in various shades of blue, sweeping across the lower half of the page. It is set against a white circular background that is itself on a light blue gradient.

ALTIUSTM
family of central venous catheters

Contents

The Altius™ Range	3
Altius™ Comprehensive Sets	4
Altius™ Basic Sets	8
Altius™ Paediatric Sets	12
Altius™ HP Comprehensive Sets	16
Guidewires	20
Insertion Packs	22

The ALTIOUS™ Range

	Multi-Tube Technology	LuerSafe™ Connector	Colour Coded Extensions and Lumens	Suture Wing	High Pressure Line
					
					

Comprehensive Sets

(Supplied in cases of 10)



ALTIUS Single Lumen Central Venous Catheter 5FR (16Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Nitinol Guidewire 0.035" x 70cm	Hydrophilic Dilator 6FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CS115501	15cm	16Ga	•	•	•	•		•	•	•	
K2CS120501	20cm		•	•	•	•		•	•	•	

ALTIUS Single Lumen Central Venous Catheter 7FR (14Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Nitinol Guidewire 0.035" x 70cm	Hydrophilic Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CS115701	15cm	14Ga	•	•	•	•		•	•	•	
K2CS120701	20cm		•	•	•	•		•	•	•	
K2CS130701	30cm		•	•	•		•	•	•	•	

ALTIUS Double Lumen Central Venous Catheter 7FR

K2CS115702	15cm	D. 16Ga P. 16Ga	•	•	•	•		•	•	•	
K2CS120702	20cm		•	•	•	•		•	•	•	
K2CS130702	30cm		•	•	•		•	•	•	•	
K2CS120702ECG	20cm		•	•	•	•		•	•	•	•
K2CS130702ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Double Lumen Central Venous Catheter 8FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Nitinol Guidewire 0.035" x 70cm	Hydrophilic Dilator 9FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CS115802	15cm	D. 14Ga P. 14Ga	•	•	•	•		•	•	•	
K2CS120802	20cm		•	•	•	•		•	•	•	

ALTIUS Triple Lumen Central Venous Catheter 7FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Nitinol Guidewire 0.035" x 70cm	Hydrophilic Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CS111703	11cm	D. 14Ga M. 18Ga P. 18Ga	•	•	•	•		•	•	•	
K2CS115703	15cm		•	•	•	•		•	•	•	
K2CS120703	20cm		•	•	•	•		•	•	•	
K2CS130703	30cm		•	•	•		•	•	•	•	
K2CS120703ECG	20cm		•	•	•	•		•	•	•	•
K2CS130703ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Triple Lumen Central Venous Catheter 8.5FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Nitinol Guidewire 0.035" x 70cm	Hydrophilic Dilator 10FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CS115853	15cm	D. 12Ga M. 14Ga P. 16Ga	•	•	•	•		•	•	•	
K2CS120853	20cm		•	•	•	•		•	•	•	

ALTIUS Quad Lumen Central Venous Catheter 8.5FR

K2CS111854	11cm	D. 14Ga M1. 16Ga M2. 16Ga P. 14Ga	•	•	•	•		•	•	•	
K2CS115854	15cm		•	•	•	•		•	•	•	
K2CS120854	20cm		•	•	•	•		•	•	•	
K2CS130854	30cm		•	•	•		•	•	•	•	
K2CS120854ECG	20cm		•	•	•	•		•	•	•	•
K2CS130854ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Quin Lumen Central Venous Catheter 8.5FR

K2CS115855	15cm	D. 14Ga M1. 18Ga M2. 18Ga M3. 18Ga P. 14Ga	•	•	•	•		•	•	•	
K2CS120855	20cm		•	•	•	•		•	•	•	
K2CS130855	30cm		•	•	•		•	•	•	•	
K2CS120855ECG	20cm		•	•	•	•		•	•	•	•
K2CS130855ECG	30cm		•	•	•		•	•	•	•	•

Basic Sets

(Supplied in cases of 10)



ALTIUS Single Lumen Central Venous Catheter 5FR (16Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Luer Slip Syringe	Altius Guidewire Dispenser	Stainless Steel Guidewire 0.035" x 50cm	Stainless Steel Guidewire 0.035" x 70cm	Dilator 6FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CV115501	15cm	16Ga	•	•	•	•		•	•	•	
K2CV120501	20cm		•	•	•	•		•	•	•	

ALTIUS Single Lumen Central Venous Catheter 7FR (14Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Luer Slip Syringe	Altius Guidewire Dispenser	Stainless Steel Guidewire 0.035" x 50cm	Stainless Steel Guidewire 0.035" x 70cm	Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CV115701	15cm	14Ga	•	•	•	•		•	•	•	
K2CV120701	20cm		•	•	•	•		•	•	•	
K2CV130701	30cm		•	•	•		•	•	•	•	

ALTIUS Double Lumen Central Venous Catheter 7FR

K2CV115702	15cm	D. 16Ga P. 16Ga	•	•	•	•		•	•	•	
K2CV120702	20cm		•	•	•	•		•	•	•	
K2CV130702	30cm		•	•	•		•	•	•	•	
K2CV120702ECG	20cm		•	•	•	•		•	•	•	•
K2CV130702ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Double Lumen Central Venous Catheter 8FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Luer Slip Syringe	Altius Guidewire Dispenser	Stainless Steel Guidewire 0.035" x 50cm	Stainless Steel Guidewire 0.035" x 70cm	Dilator 9FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CV115802	15cm	D. 14Ga P. 14Ga	•	•	•	•		•	•	•	
K2CV120802	20cm		•	•	•	•		•	•	•	

ALTIUS Triple Lumen Central Venous Catheter 7FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Luer Slip Syringe	Altius Guidewire Dispenser	Stainless Steel Guidewire 0.035" x 50cm	Stainless Steel Guidewire 0.035" x 70cm	Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CV111703	11cm	D. 14Ga M. 18Ga P. 18Ga	•	•	•	•		•	•	•	
K2CV115703	15cm		•	•	•	•		•	•	•	
K2CV120703	20cm		•	•	•	•		•	•	•	
K2CV130703	30cm		•	•	•		•	•	•	•	
K2CV120703ECG	20cm		•	•	•	•		•	•	•	•
K2CV130703ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Triple Lumen Central Venous Catheter 8.5FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Luer Slip Syringe	Altius Guidewire Dispenser	Stainless Steel Guidewire 0.035" x 50cm	Stainless Steel Guidewire 0.035" x 70cm	Dilator 10FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device	ECG Cable
K2CV115853	15cm	D. 12Ga M. 14Ga P. 16Ga	•	•	•	•		•	•	•	
K2CV120853	20cm		•	•	•	•		•	•	•	

ALTIUS Quad Lumen Central Venous Catheter 8.5FR

K2CV111854	11cm	D. 14Ga M1. 16Ga M2. 16Ga P. 14Ga	•	•	•	•		•	•	•	
K2CV115854	15cm		•	•	•	•		•	•	•	
K2CV120854	20cm		•	•	•	•		•	•	•	
K2CV130854	30cm		•	•	•		•	•	•	•	
K2CV120854ECG	20cm		•	•	•	•		•	•	•	•
K2CV130854ECG	30cm		•	•	•		•	•	•	•	•

ALTIUS Quin Lumen Central Venous Catheter 8.5FR

K2CV115855	15cm	D. 14Ga M1. 18Ga M2. 18Ga M3. 18Ga P. 14Ga	•	•	•	•		•	•	•	
K2CV120855	20cm		•	•	•	•		•	•	•	
K2CV130855	30cm		•	•	•		•	•	•	•	
K2CV120855ECG	20cm		•	•	•	•		•	•	•	•
K2CV130855ECG	30cm		•	•	•		•	•	•	•	•

Paediatric Sets

(Supplied in cases of 5)



ALTIUS Single Lumen Central Venous Catheter 2FR (22Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 4FR 5cm	Hydrophilic Dilator 4FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105201	5cm	22Ga	•	•	•	•			•		•	
K2CS108201	8cm		•	•	•	•			•		•	

ALTIUS Single Lumen Central Venous Catheter 3FR (20Ga)

K2CS105301	5cm	20Ga	•	•	•		•		•		•	
K2CS108301	8cm		•	•	•		•		•		•	
K2CS113301	13cm		•	•	•		•		•		•	

ALTIUS Single Lumen Central Venous Catheter 4FR (18Ga)

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 5FR 5cm	Hydrophilic Dilator 5FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105401	5cm	18Ga	•	•	•		•			•	•	•
K2CS108401	8cm		•	•	•		•			•	•	•
K2CS113401	13cm		•	•	•		•			•	•	•

ALTIUS Double Lumen Central Venous Catheter 3FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 4FR 5cm	Hydrophilic Dilator 4FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS104302	4cm	D. 21Ga P. 23Ga	•	•	•		•		•		•	
K2CS105302	5cm		•	•	•		•		•		•	
K2CS108302	8cm		•	•	•		•		•		•	
K2CS113302	13cm		•	•	•		•		•		•	

ALTIUS Double Lumen Central Venous Catheter 4FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 5FR 5cm	Hydrophilic Dilator 5FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105402	5cm	D. 18Ga P. 21Ga	•	•	•		•			•	•	•
K2CS108402	8cm		•	•	•		•			•	•	•
K2CS113402	13cm		•	•	•		•			•	•	•

ALTIUS Double Lumen Central Venous Catheter 5FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 6FR 5cm	Hydrophilic Dilator 6FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105502	5cm	D. 18Ga P. 20Ga	•	•	•		•			•	•	•
K2CS108502	8cm		•	•	•		•			•	•	•
K2CS113502	13cm		•	•	•		•			•	•	•

ALTIUS Triple Lumen Central Venous Catheter 4FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 5FR 5cm	Hydrophilic Dilator 5FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105403	5cm	D. 18Ga M. 23Ga P. 21Ga	•	•	•		•			•	•	•
K2CS108403	8cm		•	•	•		•			•	•	•
K2CS113403	13cm		•	•	•		•			•	•	•
K2CS130403	30cm		•	•	•	70cm				•	•	•

ALTIUS Triple Lumen Central Venous Catheter 5FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 20Ga 4cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 6FR 5cm	Hydrophilic Dilator 6FR 5cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS105503	5cm	D. 18Ga M. 22Ga P. 21Ga	•	•	•		•			•	•	•
K2CS108503	8cm		•	•	•		•			•	•	•
K2CS113503	13cm		•	•	•		•			•	•	•
K2CS130503	30cm		•	•	•	70cm				•	•	•

ALTIUS Quad Lumen Central Venous Catheter 6FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 7FR 5cm	Hydrophilic Dilator 7FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS111604	11cm	D. 16Ga M1. 20Ga	•	•	•			•		•	•	•
K2CS115604	15cm	M2. 20Ga P. 20Ga	•	•	•			•		•	•	•

ALTIUS Quin Lumen Central Venous Catheter 7FR

Contents Code	Length	Lumen Size(s)	Introducer Needle 18Ga 7cm	3ml Luer Slip Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.018" x 50cm	Nitinol Guidewire 0.021" x 50cm	Nitinol Guidewire 0.032" x 50cm	Dilator 8FR 12cm	Hydrophilic Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115705*	15cm	D. 16Ga M1. 18Ga M2. 20Ga M3. 20Ga P. 18Ga	•		•			•		•	•	•

* Please note this set includes the 5ml Guidewire Introduction Syringe instead of the 3ml Luer Slip Syringe

Comprehensive Sets

(Supplied in cases of 10)



ALTIUS HP Single Lumen Central Venous Catheter 5FR (16Ga)

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 6FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115501P	15cm	•	16Ga	•	•	•	•	•	•	•
K2CS120501P	20cm	•		•	•	•	•	•	•	•

ALTIUS HP Single Lumen Central Venous Catheter 7FR (14Ga)

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115701P	15cm	•	14Ga	•	•	•	•	•	•	•
K2CS120701P	20cm	•		•	•	•	•	•	•	•

ALTIUS HP Double Lumen Central Venous Catheter 7FR

K2CS115702P	15cm	•	D. 14Ga P. 17Ga	•	•	•	•	•	•	•
K2CS120702P	20cm	•		•	•	•	•	•	•	•

ALTIUS HP Double Lumen Central Venous Catheter 8FR

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 9FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115802P	15cm	•	D. 14Ga P. 16Ga	•	•	•	•	•	•	•
K2CS120802P	20cm	•		•	•	•	•	•	•	•

ALTIUS HP Triple Lumen Central Venous Catheter 7FR

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 8FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115703P	15cm	•	D. 16Ga M. 18Ga	•	•	•	•	•	•	•
K2CS120703P	20cm	•	P. 18Ga	•	•	•	•	•	•	•

ALTIUS HP Triple Lumen Central Venous Catheter 8FR

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 9FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115803P	15cm	•	D. 14Ga M. 18Ga	•	•	•	•	•	•	•
K2CS120803P	20cm	•	P. 18Ga	•	•	•	•	•	•	•

ALTIUS HP Quad Lumen Central Venous Catheter 8FR

K2CS115804P	15cm	•	D. 14Ga M1. 18Ga M2. 18Ga	•	•	•	•	•	•	•
K2CS120804P	20cm	•	P. 18Ga	•	•	•	•	•	•	•

ALTIUS HP Quin Lumen Central Venous Catheter 9FR

Contents Code	Length	HP Lumen (Purple)	Lumen Size(s)	Introducer Needle 18Ga 7cm	5ml Guidewire Introduction Syringe	Altius Guidewire Dispenser	Nitinol Guidewire 0.035" x 50cm	Hydrophilic Dilator 10FR 12cm	No 11 Scalpel Short Handle	Secondary Fixation Device
K2CS115905P	15cm	•	D. 16Ga M1. 20Ga M2. 20Ga	•	•	•	•	•	•	•
K2CS120905P	20cm	•	M3. 20Ga P. 20Ga	•	•	•	•	•	•	•

Guidewires

(Supplied in cases of 25)



ALTIUS Central Venous Catheter Guidewires pre-mounted on the ALTIUS Guidewire Introducer/Dispenser

Contents Code	Contents	Where used
KGW01850N	Nitinol Guidewire 0.018" 50cm	Paediatric CVCs requiring an 0.018" guidewire
KGW02150N	Nitinol Guidewire 0.021" 50cm	Paediatric CVCs requiring an 0.021" guidewire
KGW03550SW	Stainless Steel Guidewire 0.035" 50cm	Adult CVCs (inc Altius) requiring an 0.035" 50cm stainless steel guidewire
KGW03250SW	Stainless Steel Guidewire 0.032" 50cm	Adult CVCs (eg Arrow) requiring an 0.032" 50cm stainless steel guidewire
KGW03550N	Nitinol Guidewire 0.035" 50cm	Adult CVCs (inc Altius) requiring an 0.035" 50cm nitinol guidewire
KGW03250N	Nitinol Guidewire 0.032" 50cm	Adult CVCs (eg Arrow) requiring an 0.032" 50cm nitinol guidewire
KGW03570SW	Stainless Steel Guidewire 0.035" 70cm	Adult CVCs (inc Altius) requiring an 0.035" 70cm stainless steel guidewire
KGW03270SW	Stainless Steel Guidewire 0.032" 70cm	Adult CVCs (eg Arrow) requiring an 0.032" 70cm stainless steel guidewire
KGW03570N	Nitinol Guidewire 0.035" 70cm	Adult CVCs (inc Altius) requiring an 0.035" 70cm nitinol guidewire
KGW03270N	Nitinol Guidewire 0.032" 70cm	Adult CVCs (eg Arrow) requiring an 0.032" 70cm nitinol guidewire

A damaged or dropped guidewire can lead to delays in completing a procedure. Kimal offers you a full range of spare guidewires, pre-mounted in the unique Altius™ Guidewire Dispenser/Introducer for all adult and paediatric central venous catheters. The range includes a choice of size, length and either nitinol or stainless steel for adults. For paediatrics, where kinking of guidewires is an issue, we only offer nitinol.

Insertion Packs



Description		ALTIUS Basic CVC Insertion Pack	ALTIUS Intermediate CVC Insertion Pack	ALTIUS Comprehensive CVC Insertion Pack*
Contents	Quantity	K40174/CVC (Case of 25)	K40175/CVC (Case of 15)	K40176/CVC/P1 (Case of 15)
Tray 24.5cm x 13.5cm x 2.5cm	1	•	•	•
Scalpel No 11 Short Handle	1	•	•	
Scalpel No 11 (Retractable Blade)	1			•
White Swabs 10cm x 10cm	Tied 10	•	•	
Green Swabs 10cm x 10cm	Tied 5			•
120ml Gallipots	2	•	•	
Spongesticks 5cm	2	•	•	
Chloraprep® 3ml Tinted Orange**	1			•
BD 10ml Luer Slip Syringe	1	•	•	•
BD 5ml Luer Slip Syringe	1			•
Hypodermic Needle 21Ga x 3.8cm (Green)	1	•	•	•
Hypodermic Needle 23Ga x 1.6cm (Orange)	1			•
Blunt Drawing Up Needle	1			•
Mersilk Suture (Ethicon W775H)	1	•	•	•
Clear CVC Insertion Drape 65cm x 60cm with 10cm Adhesive Edge Fenestration	1	•		
Hybrid CVC Drape (Blue) 110cm x 110cm with 60cm x 60cm Clear Window with 10cm Adhesive Edge Fenestration	1		•	•
Wrapped in Tray Wrap White 75cm x 75cm	1		•	•
Hand Towels 57cm x 33cm	2		•	•
Gown SMS Large Blue	1		•	•
Wrapped in Trolley Cover Blue 90cm x 75cm	1	•		
Wrapped in Trolley Cover Blue 140cm x 140cm	1		•	•

* This pack is only for sale in the UK and Middle East

** Chloraprep™ is a registered trademark of CareFusion®

Kimal offer a choice of standard and customised insertion packs. If you have more specific requirements, please contact your Kimal Vascular Representative who will be pleased to assist you with the contents and design of your procedure pack.



**For more information, pricing or
samples speak to your local Kimal
vascular representative or contact:**

UK

Tel: +44 (0)845 437 9542

Fax: +44 (0)845 437 9544

Email: customercare@kimal.co.uk

Germany

Tel: +49 (0)2133 285 141

Fax: +49 (0)2133 285 1420

Email: info@kimal-med.de

Middle East

Tel: +97 16 552 8146

Fax: +97 16 552 8147

Email: mea@kimal.co.uk

International

Tel: +44 (0)845 437 9543

Fax: +44 (0)845 437 9544

Email: customercare@kimal.co.uk

www.kimal.com

© 2013 Kimal plc. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal plc.

Kimal plc reserve the right to change the specification of its products at any time.
Kimal plc, Arundel Road, Uxbridge, UB8 2SA.

L-0155-0062-A1