



Suitable for all Central Venous Catheters

DON'T DELAY COMPLETING A PROCEDURE...

Altius® guidewire dispensers provide kink-resistant and easy insertion of the pre-mounted guidewire.

A dropped or damaged guidewire during a CVC insertion procedure can result in delay completing the procedure, be annoying and if you do not have a spare, costly!



Insertion depth confirmed by marks on the wire and body of the dispenser


Wheel feature provides smooth controlled insertion




Kimal has a simple solution to dropped or damaged guidewires...

Altius® replacement guidewires.

We provide high quality replacement guidewires for adult patients in both stainless steel and in kink-resistant nitinol pre-mounted in our patented, unique and easy to use Altius® dispenser.



Wheel feature ensures that retracting and advancing the 'J' Tip guidewire has never been easier or smoother even with gel or blood on your gloves.



Product ordering details

Product Code	Description		
	Material	Gauge	Length
KGW03250SW	Stainless Steel	0.032"	50cm
KGW03270SW	Stainless Steel	0.032"	70cm
KGW03250N	Nitinol	0.032"	50cm
KGW03270N	Nitinol	0.032"	70cm
KGW03550SW	Stainless Steel	0.035"	50cm
KGW03570SW	Stainless Steel	0.035"	70cm
KGW03550N	Nitinol	0.035"	50cm
KGW03570N	Nitinol	0.035"	70cm

Pre-packed sterile in cases of 25

We also offer guidewires suitable for Paediatrics

We are pleased to offer replacement guidewires for use in CVC placement for paediatric patients, where we know kinking of the guidewire is seen much more. We offer nitinol guidewires pre-mounted in the dispenser in both 0.018" and 0.021" sizes.

Product Code	Description		
	Material	Gauge	Length
KGW01850N	Nitinol	0.018"	50cm
KGW02150N	Nitinol	0.021"	50cm

KIMAL
delivering healthcare innovation

Kimal
Unit 401, Pointon Way, Stonebridge Cross Business Park,
Droitwich Spa, Worcestershire, WR9 0LW
Tel: +44 (0)845 437 9542 Fax: +44 (0)845 437 9544
customer@kimal.co.uk

Kimal GmbH
Hamburger Straße 26
41540 Dormagen
Germany

Kimal International Healthcare (FZE)
PO Box 8532, SAIF - Zone P6-062, Sharjah
United Arab Emirates
Tel: + 9716 552 8146 Fax: + 971 6 552 8147
mea@kimal.co.uk

©2012 Kimal plc. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal plc. Specification: Kimal plc reserves the right to change the specification of its products/services at any time. Registered address: Kimal plc, Arundel Road, Uxbridge, UB8 2SA.

L-0155-0057-A2