



Comfy Chair C2

For Safety and Comfort during sample collection



Modern day health care providers would like to ensure comfort and safety to their patients in all stages of diagnostics and treatment. The Comfy Chair allows the patient to be seated comfortably during sample collection and also positions the arm of the patient for optimal venipuncture to the user friendliness of the technician.

Comfy Chair C2 through electrically operated synchronized movement of the back rest and leg rest allows infinite positioning of the patient to be relaxed and appropriately supported during the draw.

Should the patient feel faint during collection, you can actuate the chair into the Trendelenburg position to stabilize the patient. Movement actuation by means of wired remote.

Specifications

Seat Height	: 520mm
Width	: 550mm
Back rest length	: 820mm
Seat Length	: 545mm
Leg Rest Length	: 500mm
Profile Length	: 2000mm
Upholstery Thickness	: 2 inches
Arm Rest dimension	: 450 X 150 mm
Arm Rest Adjustment	: Height Adjustment, Lift up for easy entry of patient, angle adjustment for optimal venipuncture and may be set at an angle from the couch.
Trendlenberg position	: -10 degree.
Voltage	: 220 $\pm$ 10%
Head Rest thickness	: 4 inches.

Colour Options



QUALIMED
Assuring Quality

Plot No. 5, I D A, Veli, Trivandrum - 695 021, Kerala. India. Phone : +91 471 2501449 / 2503449
E - mail : info@qualimedindia.com Website : www.qualimedindia.com

*Specifications are likely to change due to continuous improvements