

Comfy Couch - D4

For Dialysis, Transfusion and Oncology.



QUALIMED
Assuring Quality

Functionality,
Safety and Comfort
combined with
an Elegant Design

Modern approach to healthcare demands special functional solutions for comfortable therapy of patients or donors who prefer to be treated in a position of their choice and to suit the clinical requirements of the procedure. This becomes all the more important when the patient or donor needs to be on the couch for more time. Users can choose among the different combinations provided by QualiMed suited for their procedure and price.

Benefits have been recorded in various studies of the need to eliminate the typical sickroom feeling to the patient undergoing the therapy, which facilitates faster recovery. This is possible to a large extent by

replacing the typical hospital bed, which allows very little mobility during the procedure with an attractively designed couch, which provides

various degrees of freedom to the patient or donor and the clinical staff.

Rehabilitating chronic dialysis patients are more tolerable to their repetitive and monotonous treatment if they are offered avenues for communication and other activities like carrying on a conversation, playing games, reading, writing etc. Being on a horizontal bed brings in a pessimistic attitude due to the typical sick room fixture and places limitations on the activities of the patient or donor.

The Comfy Couch is designed to provide maximum comfort and safety to the donor or patient, but also allows nurses and physicians to react in critical situations. The powerful motor system allows the shifting to shock position quickly and initiation of immediate emergency measures.



Other Features :



The Comfy Couch occupies 30 percent less space than a conventional hospital bed and hence can accommodate more number of couches in a given space.

The arm rest can be optimally set at the desired height and also can be moved to the side. This allows vascular access at a relaxed posture.

Fewer episodes of hypotension and muscle cramps due to the flexible positions which may be adapted by the patient or donor.

The height adjustment is done vertically with the help of reliable heavy duty lifting columns, thereby not causing any change in the profile of the donor.

Provided with extra thick upholstery for greater comfort during the lengthy periods of time that the treatment requires.





TECHNICAL DATA

Models	D40 – Four Motor Couch with Height Adjustment. D41 – Three Motor Couch without Height Adjustment.
Profile Length	2010mm
Seat Width	680mm
Height	600 ~ 800mm
Back Rest Adjustment	0°- 54°
Leg Rest Adjustment	0°- 60°
Trendelenburg position	-15°
Operation Method	Logically placed buttons on easy to operate remote allows quick adjustments to back rest, seat, leg rest and height. Single button emergency actuation for quick shift to shock position.
Castor Wheels	Four 100mm lockable wheels
Maximum Load	150kg
Arm Rest	Length 600mm and width 300mm Can be adjusted vertically and swivelled away and towards the body profile.
Upholstery	100mm thick upholstery with adjustable Head rest
Operating Voltage	230 ±10% AC
Tray Stand (Optional for Blood Transfusion)	Tray Stand on Castor wheels with BCM tray, Sphygmomanometer tray, Accessories tray and adjustable IV Pole.
IV Stand	Option I : Telescopic IV stand attached to side of couch and may be interchanged to either side of Couch. Option II : Telescopic IV stand on castor wheels which can be moved to either side of Couch.
Foot Support (Optional)	SS 304 with upholstery padding which can be adjusted to height of patient.
Over Bed Table(Optional)	Single Pole ABS surface on castor wheels with height adjustment (range 75 cms to 105 cms)
Handle(Optional)	SS 304 handle provided at back of couch to facilitate easy manoeuvrability
Safety Standards	Meets requirement of EN 60601-1 and EN 60601-1-2

* Specifications are likely to change due to continuous improvements

QUALIMED
Assuring Quality

Qualimed Systems, Plot # 5, Industrial Development Area, Kochu Veli, Thiruvananthapuram - 695 021,
Kerala, S. India. Tel: +91 471 2501449, Fax: +91 471 2503449, Email: sales@qualimedindia.com

www.qualimedindia.com

Colour Options



Steel Grey1001



Beige1002



Navy Blue1003



Maroon 1004